



prisma VENT50-C

Ventilator

prisma VENT50-C — The best care always.

The prisma VENT50-C ventilator is equipped with High-Flow Therapy (HFT) mode in addition to all the other features available in prisma VENT50.

In HFT mode the device provides a constant flow of five to 60 liters per minute. An oxygen feed also can be added. The flow of air is heated and humidified (for example, with prisma VENT AQUA) and delivered to the patient via a nasal cannula. The patient interface comes in different sizes for non-invasive use. A different interface is available for invasive application.

The prisma VENT50-C also may be used with a single tube with patient valve and with other tube systems with nasal cannula for High-Flow Therapy.

Like prisma VENT50, this ventilator has a pressure range up to 50 hPa. The prisma VENT50-C offers both pressure-controlled and volume-controlled mouthpiece ventilation (MPVp and MPVv). Features and integrated functions include LIAM (Lung Insufflation Assist Maneuver) for cough support and secretion management, an extensive alarm management system, rechargeable battery with a capacity of up to 12 hours and an integrated oxygen port.

High-Flow Therapy (HFT) mode

- Washes out nasopharyngeal dead space and thereby eliminates CO₂
- Improves mucociliary clearance by humidifying and warming the upper airways
- Improves oxygenation
- · Reduces spontaneous breathing frequency

Brilliant display

• Large display, high contrast, clear and natural colors

Simple operation and setting process

- Intuitive menus, simple use, fast optimization
- Simple setting with 360° easy turn knob and quick selection key

Highly compatible

- · Suitable for a range of tube systems
- Tube system for High-Flow Therapy
- Supplies power to additional devices
- · Works with commercially available cable



Article Nr.

prisma VENT50-C

WM 30680-1110

prisma VENT50-C with integrated rechargeable battery

WM 30690-1110

Please refer to the currently valid price list.

Accessories

prisma VENT AQUA 100506

Respiratory gas humidifier

For High-Flow Therapy only

271708 Single tube system with patient valve

High-Flow leakage

271705

heated (i), with auto-fill chamber for prisma VENT AQUA

(i), auto-fill chamber, passive valve, Connection for HFT nasal cannula (150 cm + 60 cm, 22 mm Ø)

(150 cm + 60 cm, 22 mm Ø)

Oxi Plus High-Flow nasal cannula

Single tube system with patient valve heated (i), with auto-fill chamber for prisma VENT AQUA (150 cm + 60 cm, 15 mm Ø)

Sold in packages of 15 Size 2 270822 270823 Size 3 270824 Size 4

You'll find the complete line of accessories in the prisma VENT accessories catalog. The prisma VENT50-C is compatible with existing prisma VENT accessories.

271707

Tochnical data

| Technical data | |
|---|--|
| Product Class as per 93/42/EEC | lla |
| Dimensions $W \times H \times D$ in cm | 21.8 × 17.5 × 21.8 |
| Weight | 2.5 kg |
| Weight of internal rechargeable battery | 0.63 kg |
| Maximum air flow at 20 hPa | >220 l/min |
| Operating time of internal rechargeable battery with the following settings: T mode, $f = 20$ /min, $Ti = 1$ s, PEEP = 4 hPA, $Vt = 800$ ml Passive lungs: Resistance R = 5 hPa (I/s); Compliance C = 50 ml/hPa | > 10 hours |
| IPAP pressure range | 4 to 50 hPa |
| PEEP pressure range | 4 to 25 hPa (leakage system) 0 to 25 hPa (valve system) |
| Respiratory rate | 0 to 60 1/min |
| Target volume | 100 ml to 2000 ml |
| Trigger levels | |
| Inspiration | I (high sensitivity) to 8 (low sensitivity) |
| Expiration | 95% to 5% of maximum flow in steps of 5% |
| Pressure increase speed | Level 1: 100 hPa/s Level 3: 50 hPa/s Level 2: 80 hPa/s Level 4: 20 hPa/s |
| Pressure decrease speed | Level 1: 100 hPa/s Level 3: 50 hPa/s Level 2: 80 hPa/s Level 4: 20 hPa/s |
| Tidal volume | 100 ml to 2000 ml |
| Maximum allowed flow with oxygen feed | 15 l/min |

Löwenstein Medical Technology Kronsaalsweg 40, 22525 Hamburg

info@loewensteinmedical.de www.loewensteinmedical.de

Germany T: +49 40 54702-0 F: +49 40 54702-461



medical

© Copyright protected. Reproduction of any type only with express permission from Löwenstein Medical Technology. We are not responsible for printing or typographical errors.