

PLUS A7

ELECTROMECHANIC ICU BED

 **Medikal 2000**

Technical Specifications

- Electronic adjustment of backrest, leg rest, height-up & down, Fowler positions
- remote control
- Rastomat Mechanism
- 4- piece, HPL lying surface
- Headboard & footboard removable, side rails are foldable
- Plastic bumpers in corners
- Drainage bag holder hooks.
- Holes for accessories in all corners
- Height Adjustable IV Pole
- 11 cm thick mattress compatible with bed lying surface
- X-ray transparent lying surface
- Nano ceramic, electrostatic epoxy polyester powdermetal parts
- Backrest regression
- Minimum IPX4 protection
- Angle indicators on the side rails
- 3 linear motors



— Position 1 —



— Position 2 —



— Position 3 —



— Position 4 —



— Position 5 —

Options

- Trendelenburg Or Reverse Trendelenburg by Gas spring
- Back up battery for non-stop operation in case of power failure
- Manual and electronically CPR (Reset).
- Nurse Control Panel
- Colored sticker options for head / foot boards and side rails
- X-Ray cassette carrier.
- Side Rail Control Panels