







# PLUS A8

ELECTROMECHANIC  
ICU BED

 Medikal 2000



# PLUS A8

## ELECTROMECHANIC ICU BED



### Technical Specifications

- Electronic adjustment of backrest, leg rest, height-up & down, Trendelenburg, Reverse Trendelenburg, Cardiac and Fowler positions
- Side Rail Control Panels
- Rastomat Mechanism
- 4- piece, HPL lying surface
- Headboard & footboard removable, side rails are foldable
- 360° turning swivel casters equipped with central and directional lock
- Plastic bumpers in corners
- Drainage bag holder hooks.
- Holes for accessories in all corners
- Height Adjustable IV Pole
- 11 cm thick mattress compatible with bed lying surface
- X-ray transparent lying surface
- Nano ceramic, electrostatic epoxy polyester powder metal parts
- Backrest regression
- Minimum IPX4 protection
- Angle indicators on the side rails
- 4 linear motors



### Options

- Back up battery for non-stop operation in case of power failure
- Manual and electronically CPR (Reset).
- Nurse Control Panel
- Colored sticker options for head / foot boards and side rails
- X-Ray cassette carrier.