






M2-T800

HYDRAULIC
EMERGENCY STRETCHER



M2-T800

HYDRAULIC EMERGENCY STRETCHER



Specifications :

- 230 kg load capacity.
- hydraulic controlled height, Trendelenburg and reverse Trendelenburg adjustments with foot pedals
- hydraulic system designed with aluminum columns.
- back rest adjustment controlled with gas springs.
- 200 mm diameter castors capable of turning 360°.
- the fifth wheel support for 360-degree turn in tight spaces.
- Central break system with corridor system
- brake pedal located on the four corners.
- compact laminate lying surface.
- easily cleanable, two-piece compact laminate lying surface.
- a full-length x-ray permeability.
- x-ray cassette holder.
- antibacterial seamless mattress.
- easy foldable side rails.
- Metal parts with zinc phosphate and varnished with electrostatic epoxy polyester powder.
- IV pole inserts and protective bumpers on the four corners.
- push bar manufactured from stainless steel.
- ABS plastic covered lower chassis.
- designed ABS plastic to carry oxygen tube on.

