



Nano

Portable Oxygen Concentrator



Model Part Number

855 100 to 240 VAC
50 to 60Hz

Specifications

Average Oxygen Content	90% (+6% / -3%) at all settings
Rechargeable Battery	14.4 VDC Voltage / 6.8Ah Rated Capacity
Breathing Frequency	10 to 40 BPM
Weight	4.7 lbs / 2.1 kg
Dimensions	8.3"W x 3.5"D x 6.5"H 23 cm x 9 cm x 17 cm

2010-8100 REV A

Light Weight!

Take it Anywhere!

Bright LCD Screen for Easy Viewing!

Easy Touch Button Controls!

Convenient Carrying Bag Included!

NIDEK

Medical



MADE IN
USA

CE
0413